

## Supplementary Files

**Figure S1** presents the representative isotype control staining (gray) for each marker used to detect developmental marker expression in neonatal (blue) and adult (red) CPCs.

**Figure S2** outlines the directed differentiation assays for cardiomyocyte, endothelial, and smooth muscle pathway induction (A). Cardiomyocyte (B-C), endothelial cell (D-E), endothelial tube (F), and smooth muscle cell (G-H) induction was screened after 14 days by RT-PCR and flow cytometry. n=3-9 replicates. Data: mean  $\pm$  SEM, \*p<0.05, \*\*p<0.01, \*\*\*p<0.001.

**Figure S3** provides a schematic of the experiment, including how cell culture systems were prepared prior to launch and a timeline of experiments once aboard the ISS.

**Figure S4** indicates changes in cardiac (neonatal, A; adult, B) and oxidative (neonatal, C; adult, D) stress response as well as apoptosis activation (neonatal, E; adult, F) in adult and neonatal CPCs flown aboard the ISS for 12 days. n=3 biological replicates. Data: mean  $\pm$  SEM, \*p<0.05, \*\*p<0.01, \*\*\*p<0.001.


**Figure S5** shows the cell count (A) and viability (B) of adult and neonatal CPCs upon return to Earth after 30 days of culture aboard the ISS. n=4 pooled clones measured in technical replicate for cell count and viability results. Data are reported as the mean  $\pm$  SEM.

**Table S1** reports KEGG analysis of targets of significantly altered microRNAs by age.

**Table S2** indicates the manufacturer, isotype, concentration, species, clone, catalog number, and lot number of all antibodies used in all experiments.

**Table S3** details primer pair sequences (reported 5' to 3') used in RT-PCR experiments as well as the catalog number, manufacturer, and name of microRNA primers that were not originally included in microarray plates.

**Figure S1**



**Figure S2**





Figure S3


Biocell Seeding




Experiment Timeline

Day 0	Day 3	Day 8	Day 12	Day 16	Day 20
Launch (37°C & 5% CO <sub>2</sub> )	Arrives to ISS  Media Change (Groups A & B)	Media Change (Groups A & B)	Fixation (Biocell Group A)  Media Change (Biocell Group B)	Media Change (Biocell Group B)	Media Change (Biocell Group B)
Day 24	Day 27	Day 30	Day 31 – Return to Lab		
Media Change (Biocell Group B)	Media Change (Biocell Group B)	Return to Earth (37°C & 5% CO <sub>2</sub> )	<u>Group A</u> RNA Purification  <u>Group B</u> Endothelial-like Tube Assay Cell Cycle Analysis Migration Assay Protein Isolation		

Figure S4



**Figure S5**



**Table S1****Neonatal CPC KEGG Analysis Results**

<b>KEGG pathway</b>	<b>P-Value</b>	<b>Genes</b>	<b>miRNAs</b>
Fatty acid biosynthesis	0	2	1
ECM-receptor interaction	0	14	2
Hippo signaling pathway	1.94E-13	43	3
Cell cycle	1.56E-05	43	3
Lysine degradation	2.80E-05	16	2
Glycosphingolipid biosynthesis - lacto and neolacto series	0.000139519	2	2
p53 signaling pathway	0.002223088	25	3
Protein processing in endoplasmic reticulum	0.02048608	51	2
Endocytosis	0.02304169	52	2
Thyroid hormone signaling pathway	0.04192693	34	2

**Adult CPC KEGG Analysis Results**

<b>KEGG pathway</b>	<b>P-Value</b>	<b>Genes</b>	<b>miRNAs</b>
Fatty acid biosynthesis	0	4	2
ECM-receptor interaction	0	19	4
Hippo signaling pathway	0	87	8
Lysine degradation	1.11E-16	28	6
Adherens junction	1.26E-10	51	5
Cell cycle	2.97E-09	72	3
TGF-beta signaling pathway	3.00E-08	49	5
Protein processing in endoplasmic reticulum	4.16E-08	99	4
Fatty acid metabolism	3.36E-06	16	2
p53 signaling pathway	2.40E-05	41	4
Endocytosis	0.000157412	76	3
FoxO signaling pathway	0.000502794	51	3
RNA transport	0.0013397	73	2
Thyroid hormone signaling pathway	0.003808595	48	3
Estrogen signaling pathway	0.003832405	29	3
Ubiquitin mediated proteolysis	0.00920295	70	3
Signaling pathways regulating pluripotency of stem cells	0.009475487	56	3
Glycosaminoglycan biosynthesis - chondroitin sulfate / dermatan sulfate	0.01175345	9	2



**Table S2****Antibodies used for developmental marker detection**

Antibody	Manufacturer	Isotype	Conc.	Species	Clone	Cat. No.	Lot No.
Mesp1-DyLight 405	Novus Biologicals	IgG	0.8 mg/mL	Rabbit	2030B	MAB9219V	CKNQ01-1-050317-V
Isotype-DyLight405	Novus Biologicals	IgG	1.05 mg/mL	Rabbit	NBP2-36463V	NBP2-36463V	31908-050317-V
Viability 405/520	Miltenyi Biotec	n/a	n/a	n/a	n/a	130-109-814	5170315796
Islet1	Abcam	IgG1	0.83 mg/mL	Mouse	1H9	ab86472	GR273015-3
Fluorescein Labeling Kit	Novus Biologicals	n/a		n/a	10934	707-0030	1BS3384
Isotype-FITC	BioLegend	IgG1 kappa	0.5 mg/mL	Mouse	MOPC-21	400107	B199152
PDGFR $\alpha$ -PE	BioLegend	IgG1 kappa	100 $\mu$ g/mL	Mouse	16A1	323505	B192368
Isotype-PE	BioLegend	IgG1 kappa	200 $\mu$ g/mL	Mouse	MOPC-21	400113	B214532
CXCR4-PE/Vio770	Miltenyi Biotec	IgG1	82.5 $\mu$ g/mL	n/a	REA649	130-109-887	5170503056
REA Ctrl-PE/Vio770	Miltenyi Biotec	IgG1	20 $\mu$ g/mL	n/a	REA293	130-104-616	5170201559
cKit-DyLight650	Novus Biologicals	IgG2B kappa	500 $\mu$ g/mL	Mouse	2B8	NB100-77477C	B147020-A
Isotype-AlexaFluor 647	R&D Systems	IgG2B	10 $\mu$ g/mL	Rat	141945	IC013R	AEIU0114121
SSEA1-APC/Vio770	Miltenyi Biotec	IgG1	8.25 $\mu$ g/mL	n/a	REA321	130-104-992	5161207356
REA Ctrl-APC/Vio770	Miltenyi Biotec	IgG1	30 $\mu$ g/mL	n/a	REA293	130-104-618	5170201560

**Antibodies used for directed differentiation assays**

Antibody	Manufacturer	Isotype	Conc.	Species	Clone	Cat. No.	Lot No.
TropT-PE	Miltenyi Biotec	IgG1	10 $\mu$ g/mL	n/a	REA400	130-106-746	5160426814
Isotype-PE	BD Pharmingen	IgG1 kappa	0.5 mg/mL	Mouse	MOPC-21	556650	19587
vWF-AlexaFluor 657	Novus Biologicals	IgG1 kappa	0.75 mg/mL	Mouse	3E2D10	NBP2-34535AF647	7450-1PABX160519120816-AF647
Isotype-AlexaFluor 647	R&D Systems	IgG2B	10 $\mu$ g/mL	Rat	141945	IC013R	AEIU0114121
$\alpha$ SMA-PerCP	Novus Biologicals	IgG2a kappa	0.4 mg/mL	Mouse	1A4/a sm1	NBP2-34522PCP	59-1PABX170619-090617-PCP
Isotype-PerCP	Miltenyi Biotec	IgG1	55 $\mu$ g/mL	Mouse	IS5-21F5	130-094-968	5110906188

**Table S3**

Gene	Forward Sequence (5' to 3')	Reverse Sequence (5' to 3')
ACTA2	AGCTTTCAGCTTCCCTGAACA	TACAGAGCCCAGAGCCATTG
ACTIN	TTTGAATGATGAGCCTTCGTCCCC	GTCTCAAGTCAGTGTACAGGTAAGC
ATM	GGGCGAGCCGCAAACGCTAA	TTCGGCCCGTCGGAGCAAAC
E2F1	GACCATCAGTACCTGGCCGAGAG	GACGACACCGTCAGCCGAGTG
HGF	CACGAACACAGCTTTTTGCC	TGATCCCAGCGCTGACAAAT
HMOX1	CTCTCGAGCGTCCTCA	ACTATCAGACAATGTTGT
HSP70	TGACCAAGATGAAGGAGATCG	GTCAAAGATGAGCACGTTGC
IGF1	CAGAGCAGATAGAGCCTGCG	CAGGTAACCTCGTGACAGAGCA
MESP1	TAGGCCTCAGCGAGGAGAGT	TCCCTTGTCACTTGGGCTCC
MLC2V	GGTGCTGAAGGCTGATTACGTT	TATTGGAACATGGCCTCTGGAT
MYH11	CAAATACGCGGATGAGAGGGA	CTCATGGACGTTCTTGCCCA
NKX2-5	CGCCGCTCCAGTTCATAG	GGTGGAGCTGGAGAAGACAGA
PECAM1	AACGGAAGGCTCCCTTGATG	TAAGAACCGGCAGCTTAGCC
POU5F1	AACCTGGAGTTTGTGCCAGGGTTT	TGAACTTCACCTTCCCTCCAACCA
RAD23A	GTATCGGAGCAGCCGGCCAC	TCCCAGGGGGCTCGTTCAG
RAD50	CTACGGCTTTGCGTCCCCGG	ACACCAGCTTTCCCCGC
SDF1A/CXCL12	CTACAGATGCCCATGCCGAT	GTGGGTCTAGCGGAAAGTCC
SOX2	AACCAGCGCATGGACAGTTA	GACTTGACCACCGAACCCAT
TERT	AGAGTGTCTGGAGCAAGTTGC	CGTAGTCCATGTTTACAATCG
TNNT2	GTGGGAAGAGGCAGACTGAG	ATAGATGCTCTGCCACAGC
TP53INP1	CACAACAACAAAAGGACTTGGACT	TTGAGCTTCCACTCTGGGAC
TXNIP	TCAGTATTGCAGGGCTTGGC	GTCTCTTGAGTTGGCTGGCT
VEGFA	CAGCGAAAGCGACAGGGGCA	GCTGGAGCACTGTCTGCGCA
VWF	ACACCTGCATTTGCCGAAAC	ATGCGGAGGTCACCTTTCAG
YAP1	TCCCAGATGAACGTCACAGC	TCATGGCAAACGAGGGTCA

Symbol	MicroRNA	Catalog Number	Manufacturer
hsa-miR-99a-5p	Hs_miR-99a_2	MS00032158	Qiagen
hsa-miR-100-5p	Hs_miR-100	Ms00031234	Qiagen
SNORD72	Hs_SNORD72_11	MS00033719	Qiagen
SNORD96A	Hs_SNORD96A_11	MS00033733	Qiagen